

Titan OB/GYN Stretcher



Key Features:

- 550 lbs (250 kg) weight capacity
- 21.7" (55 cm) low height
- Storage area for patient belongings
- Radiolucent fowler section
- Easily accessible brake pedals
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces
- Single-hand-operated lowering leg section
- Retractable fluid basin holder

The Amico Titan OB/GYN Stretcher is designed with functionality and convenience in mind, allowing caregivers to provide a more meaningful level of care to patients.

This product offers key features such as a one hand operated lowering leg section, a deep well storage area for patient belongings and an ergonomic design to help improve efficiency of healthcare. To accommodate obstetrics patients with limited mobility, the Titan OB/GYN Stretcher is built on Amico's best-in-class, low height stretcher platform.

The Titan OB/GYN Stretcher can also serve as an interim transport stretcher for moving patients between departments.

Standard Features:

- Two IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Fold-down siderails
- 8" (20 cm) wheel casters
- Two side mounted pedals for patient surface height adjustment
- Hydraulic lift system capable of supporting up to 550 lbs (250 kg)
- Central brake locking system
- Ergonomic push handles
- Side bumpers
- Head and foot-end mattress retainer
- Retractable fluid basin holder

Model Number & Approval

Model Number	S-AM1-3005
Agency Approval	cULus

Overall Dimensions

Overall Length	84" (213 cm)
Overall Width	32.5" (82.6 cm)

Litter Positioning

Trend/Reverse Trend	15.5°-17°
Maximum Fowler Section Inclination	81.5°
Maximum Knee Section Inclination (optional)	24°

Height Range

Lowest	21.7" (55 cm)
Highest	33" (83.8 cm)

Surface Dimensions

Surface Width	25" (63.5 cm)
Surface Length	77.5" (197 cm)
Surface Thickness	4" (10 cm)

Weight Capacity & Caster Dimension

Maximum Weight Limit	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Foley bag holder
- Disposable fluid basins
- Overbed tables
- Permanent IV pole
- Stirrups (calf supports)

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: Three years parts

Amico Beds Corporation 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 71 East Industry Court, Deer Park, NY 11729, USA
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862